- DE Gebrauchsanweisung
- GB Instructions for use
- FR Mode d'emploi
- ES Instrucciones de uso
- IT Istruzioni per l'uso
- SE Bruksanvisning


Manufacturer:
pro3dure medical GmbH
Am Burgberg 13
D - 58642 Iserlohn
Tel. +49 (0)2374 920050-10
info@pro3dure.com
www.pro3dure.com

Distributor:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

Distributed in USA / Canada exclusively by:
Kulzer, LLC
4315 South Lafayette Blvd.
South Bend, IN 46614-2517
1-800-431-1785

Caution: Federal law restricts this device to sale by or on the order of a dentist.

Medizinprodukt
Medical device
Dispositif médical
Producto sanitario
Dispositivo medico
Medicinteknisk produkt

Explanation of symbols on labelling	
	Batch code
	Manufacturer
	Keep away from sunlight
	Consult instructions for use
	Use-by date
	Storage temperature

99001509/00

GB Instructions for use dima Print Splint clear

Product description

dima Print Splint Clear Guide is a resin for the 3D printing of dental objects based on image projection systems (≤ 405 nm). The formulation of dima Print Splint Clear Guide is optimized for the requirements of a robust production guaranteeing constant high quality. The dima Print Splint Clear Guide is biocompatible and additionally meets all mechanical and application demands. The material can be used for build processes with layer thicknesses from 25 up to 100 μ m. It is recommended to use a compatible light curing device (e.g. Hilite Power 3D) from Kulzer for post curing.

Processing

- dima Print Splint Clear Guide bottles should be agitated before use for at least 5 minutes.
- Make sure that dima Print Splint Clear Guide material is temperature adjusted up to 23°C to 30°C.
- Carefully pour dima Print Splint Clear Guide material to the filling line into the vat of the image projection unit.
- Bubbles can be removed with a cleaned spatula or by a recoater routine.
- Make sure you use the correct parameters for dima Print Splint Clear Guide and follow the instructions for use for cara print 4.0.
- After the build process is finished a direct post treatment is recommended. If this cannot be guaranteed, leave the produced objects in cara Print 4.0 until post-processing.
- Cleaning is carried out in suitable cleaning solution (e.g. isopropanol) in a 2-step procedure by pre-cleaning for 3 minutes and post-cleaning in a clean bath for 2 minutes, in an ultrasonic bath. After cleaning, the objects are cured in a suitable and compatible light curing device with the appropriate program for dima Print Splint Clear Guide. The light curing programs can be found on our homepage www.kulzer.com
- The dental objects generated out of the generative resin dima Print Splint Clear Guide can be coated and repaired as usual.
- Impurities due to operation mistakes cannot be excluded. With respect to the low viscosity of the resin, it is possible to filter the dima Print Splint Clear Guide. It is recommended to filter and stir up the resin on a regular basis. Further information on processing can be found in the cara Print 4.0 application guide at www.kulzer.com

Important

- To avoid detrimental effects on material quality, do not expose the liquid material to irradiation under any circumstances.
- Deviations from the described manufacturing process may lead to different mechanical and optical properties of the dima Print Splint Clear Guide material.
- Ensure personal protective gear is worn during processing.
- Caution: Once polymerised, resins are chemically resistant - avoid contact with clothing!
- Avoid any contact with skin and eyes. In case of accidental contact, rinse with adequate running water, consult a doctor if necessary.
- The lot number and the best before date are indicated on each dima Print Splint Clear Guide packaging. For inquiries, please always indicate the lot number of the product. Do not use the product after the "best before" date.

Safety advice

Neither Kulzer or pro3dure medical GmbH are liable for any damages caused by improper application of the material. To be used by trained specialist personnel for the purpose indicated only.

Technical data:

- Color: clear transparent
- Density: ca. 1.1 g/ml
- Viscosity: ca. 0,7 Pa s
- Post cured material: (depends on postcuring unit) Elastic modulus: ca. 1700 MPa; Flexural strength: ca. 90 MPa Elongation at break: ca. 10 % Hardness: ca. 80 Shore D

Storage:

According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

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