


IT'S MORE THAN
CLEANING.



NEW



IT'S TIME SAVING.

CARA PRINT CLEAN UP

Save time and effort with the cara Print CleanUp, your reliable solution for efficiently cleaning 3D-printed dental models with ease.

Designed for versatility, it works seamlessly with multiple cleaning buckets and offers hands-free operation to clean parts directly on the build platform.

Audible and light alerts ensure you stay in control without constant monitoring.

Using isopropyl alcohol (IPA), the cara Print CleanUp delivers thorough, precise cleaning with minimal manual intervention – making it the perfect addition to your workflow.



[kulzer.com/
mycleanup](https://kulzer.com/mycleanup)



KULZER
MITSUI CHEMICALS GROUP

CARA PRINT CLEAN UP

- ▶ **Versatile use:** One base with multiple cleaning buckets.
- ▶ **User Notifications:** Audible and light signals upon completion.
- ▶ **Hands-Free Operation:** Clean parts without removing them from the build platform.

Main parameters

Size (l x w x h) [mm]	255 x 155 x 385.5
Maximum washing volume (w/o build platform) (l x w x h) [mm]	201 x 124 x 255
Weight [kg]	2.27
Rated power [W]	12
Input Voltage	110V/240V 50/60Hz
Time settings [min]	0–30 in 30sec intervals



NEW

NEW