





## DIMA PRINT PERFORMANCE OVERVIEW


PICTURE	MATERIAL	INDICATION	EU MD CLASS	COLOUR	PRINTING TIME	EQUIVALENT PRINTING TIME PER PART *	COST PER PART **
	dima Print C&B temp	Temporary crowns & bridges	Ila	A1, A2, A3, A3.5, B1, B2, BL1	1 h 34 min (for 110 units)	3 min	0,38 €
	dima Print Cast ruby	CAD to cast	not applicable	red	1 h 28 min (for various CAD-to-Cast objects)	5 min	0,28 €
	dima Print Denture Base	Full and partial removable dentures and baseplates	Ila	original pink, light pink, light reddish pink, dark pink	1 h 44 min (for 10 arches)	12 min	5,91 €
	dima Print Denture Teeth	Denture Teeth	Ila	A1, A2, A3, A3.5, B1, B2, BL1	1 h 35 min (for 8 arches)	4 min	0,28 €
	dima Print Splint clear	Splint / Nightguards	Ila	transparent clear	1 h 18 min (for 6 Splints)	13 min	2,25 €
	dima Print Stone beige	Dental Models	not applicable	gypsum beige	1 h 30 min (for 4 big models)	22 min	12,47 €

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

\* depending on various factors, especially on the utilisation of the build platform; calculation based on build platform utilisation as shown in picture

\*\* based on the RRP of the respective dima Print material only

Find the current version here: [kulzer.com/mycube-support](https://kulzer.com/mycube-support)



**KULZER**  
MITSUI CHEMICALS GROUP

## DIMA PRINT PERFORMANCE OVERVIEW

PICTURE	MATERIAL	INDICATION	EU MD CLASS	COLOUR	PRINTING TIME	EQUIVALENT PRINTING TIME PER PART *	COST PER PART **
	dima Print Stone beige	Dental Models	not applicable	gypsum beige	1 h 18 min (for 5 models)	13 min	9,54 €
	dima Print Stone beige	Dental Models	not applicable	gypsum beige	1 h 15 min (for 1 model)	1 h 15 min	9,54 €
	dima Print Guide & Tray	Impression Trays	I	transparent light orange	1 h 43 min (for 5 trays)	28 min	6,25 €
	dima Print Guide & Tray	Surgical Guides	I	transparent light orange	57 min (for 6 big guides)	10 min	2,32 €
	dima Print Denture Base Try-In	Try-ins	IIa	white	1 h 01 min (for 6 arches)	10 min	9,79 €
	dima Print Stone beige Smart	Dental Models	not applicable	gypsum beige	21 min (for 1 model for splints)	21 min	3,53 €

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

\* depending on various factors, especially on the utilisation of the build platform; calculation based on build platform utilisation as shown in picture

\*\* based on the RRP of the respective dima Print material only

Find the current version here: [kulzer.com/mycube-support](https://kulzer.com/mycube-support)



## DIMA PRINT PERFORMANCE OVERVIEW

PICTURE	MATERIAL	INDICATION	EU MD CLASS	COLOUR	PRINTING TIME	EQUIVALENT PRINTING TIME PER PART *	COST PER PART **
	dima Print Stone beige Smart	Dental Models	not applicable	gypsum beige	23 min (for 6 models for splints)	4 min	3,53 €
	dima Print Denture Base 2	Full and partial removable dentures and baseplates	IIa	orange-pink	1h 54 min (for 10 arches)	13 min	5,07 €
	dima Print Gingiva Mask	Gingiva Masks	not applicable	pink	17 min (for 10 gingiva masks)	2 min	1,42 €
	dima Print Stone teal	Dental models for deep drawing	not applicable	turquoise	1h 57 min (for 5 models)	23 min	10,30 €
	dima Print Denture Teeth 2	Temporary Crowns & Bridges	IIa	A1, A2, A3, A3.5, B1, B2, BL1	1h 59 min (for 100 units)	4 min	0,93 €
	dima Print Denture Teeth 2	Denture Teeth	IIa	A1, A2, A3, A3.5, B1, B2, BL1	1h 56 min (for 8 arches of teeth)	5 min	1,12 €

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.




\* depending on various factors, especially on the utilisation of the build platform; calculation based on build platform utilisation as shown in picture

\*\* based on the RRP of the respective dima Print material only

Find the current version here: [kulzer.com/mycube-support](https://kulzer.com/mycube-support)



## DIMA PRINT PERFORMANCE OVERVIEW

PICTURE	MATERIAL	INDICATION	EU MD CLASS	COLOUR	PRINTING TIME	EQUIVALENT PRINTING TIME PER PART *	COST PER PART **
	dima Print Teeth & Temp	Temporary Crowns & Bridges	Ila	A1, A2, A3, A3.5, B1, B2, BL1	43 min (for 100 units)	1 min	0,98 €
	dima Print Teeth & Temp	Denture Teeth	Ila	A1, A2, A3, A3.5, B1, B2, BL1	38 min (for 8 arches of teeth)	1 min	1,39 €
	dima Print Stone gray	Dental Models	not applicable	gray	1 h 30 min (for 4 big models)	22 min	16,52 €

Availability in your country may be restricted due to legal requirements. Please check the availability for your country with Kulzer Customer Service Center.

\* depending on various factors, especially on the utilisation of the build platform; calculation based on build platform utilisation as shown in picture

\*\* based on RRP

Find the current version here: [kulzer.com/mycube-support](https://kulzer.com/mycube-support)

