

Dri-Clav VK-4 Ingredient List



VK-4 Tartar and Stain Remover is a ready-to-use formula especially for removal of calculus, tartar, tobacco and food stains from dentures, bridges and orthodontic appliances.

Use undiluted VK-4 in a covered glass beaker immersed in VK-1 solution that acts as the ultrasonic conductor between tank and beaker. Cover item with VK-4 chemical; clean 5-10 minutes; rinse item with water. Dentures, bridges, orthodontic appliances and other acrylic items should be neutralized in a solution of 1 teaspoon bicarbonate of soap diluted in a cup of water.

Product Name: Dri-Clav VK-4 Tartar & Stain Remover

Product Code:

Ingredient	CAS No.	Function
Water	7732-18-5	Solvent
Sulfamic Acid	5329-14-6	Acidic cleaning agent

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP