

VacuKleen E2 Ingredient List



Cleans with a lemon flash

VacuKleen E2 is an Enzymatic Evacuation System Cleaner that delivers the strength of enzymatic cleaning power for removing tough deposits.

VacuKleen E2 is effective for:
Easy & Efficient Cleaning of your Evacuation System

VacuKleen E2 is available:
For Quarts: Box of 60 Pouches or Box of 140 Pouches
For Gallons: Box of 25 Pouches

Product Name: VacuKleen

Product Code:

Ingredient	CAS No.	Function
Sodium Bicarbonate	144-55-8	Alkalizer, buffer
Sodium Carbonate	497-19-8	Alkalizer, buffer
Pluronic F68	9003-11-6	Non-ionic surfactant
Citric Acid	77-92-9	Preservative, effervescing agent, pH control agent
Deterzyme APUG 380	9014-01-1	Protease enzyme
Tetrasodium ethylenediaminetetraacetate	64-02-8	Water softener, chelating agent
Sodium Benzoate	532-32-1	Preservative
Lemon Powder Fragrance	1180-71-8	Fragrance ingredient
FD&C Yellow#5	1934-21-0	Colorant

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP