



VK Solutions and Tablets®

Powerful cleaning with minimal effort.

- Powerful ultrasonic energy quickly removes and breaks down hard to clean deposits like blood, bone, and food particles that chemical detergents can leave behind.
- Powerful, yet gentle on your instruments.
- Biodegradable and pH neutral*.
- Non-toxic, non-corrosive formula contains no chlorine or harsh alkaline to damage your instruments*.

*Refers to VibraKleen E2 Ultrasonic Cleaning Tablets.

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

Your Complete Ultrasonic Cleaning System

Product	Features	Benefits	Uses	Packaging	Product Numbers
VIBRAKLEEN E²® ULTRASONIC CLEANING SOLUTIONS					
VibraKleen E² Ultrasonic Cleaning Solutions	<ul style="list-style-type: none"> • Effervescent, enzymatic cleaning tablets • Non-corrosive and low-foaming • Biodegradable, non-toxic & pH neutral • No chlorine or harsh alkaline • Drop directly into ultrasonic cleaner • Compact, convenient package of 20 tablets dilutes to 10 gallons • Clean mint scent • Effective as evacuation cleaner • Power of protease enzymes 	<ul style="list-style-type: none"> • Breaks down hard-to-clean deposits in ultrasonic and evacuation systems • Easy on instruments • Environmentally friendly and safe to work with • Won't damage stainless steel equipment • No messy mixing of liquid concentrates • Saves on counter and storage space • Minimizes any odor from tank, pleasant to work with • One product for multiple uses • Effective, fast removal of blood, bone and food particles that chemical detergents leave behind 	<ul style="list-style-type: none"> • For general purpose ultrasonic cleaning, evacuation system cleaning and pre-soak 	Box of 20 Box of 80 2 tablets = 1 gallon	50037100 — Box of 20 (Case Lot 24) 50037101 — Clinic Pack— Box of 80 (Case Lot 12)
DRI-CLAVE® VK ULTRASONIC CLEANING SOLUTIONS					
VK-1 General Purpose Cleaner	<ul style="list-style-type: none"> • Efficient cleansing time (2–10) minutes • 10:1 concentrated formulation • Color coded solution • Anticorrosive 	<ul style="list-style-type: none"> • Protects your equipment • Avoid messy cleaning • Economical gallon container • Mistake-free identification 	<ul style="list-style-type: none"> • Removes debris from instruments, burs, general equipment, ortho appliances, gold and porcelain • Use directly in ultrasonic tank 	Gallon— blue Liquid	50036202* 4000239* (Canadian customers only)
VK-3 Plaster and Stone Remover	<ul style="list-style-type: none"> • Efficient cleansing time (5–10) minutes • No diluting necessary—use at full strength • Color coded solution 	<ul style="list-style-type: none"> • Protects your instruments • Avoid messy cleaning • Saves time • Economical gallon container • Mistake-free identification 	<ul style="list-style-type: none"> • For removal of plaster and stone • Do not use with aluminum • Place in glass breaker inside ultrasonic cleaner 	Gallon— pink Liquid	50036206* 4000240* (Canadian customers only)
VK-4 Tartar and Stain Remover	<ul style="list-style-type: none"> • Efficient cleansing time (5–10) minutes • No diluting necessary—use at full strength 	<ul style="list-style-type: none"> • Avoid messy cleaning • Ready-to-use formula • Saves time • Economical gallon container 	<ul style="list-style-type: none"> • Removal of tartar and stains • 5–10 minutes only • Place in glass beaker inside ultrasonic cleaner 	Gallon— clear Liquid	50036208* 4000241* (Canadian customers only)
VK-6 Temporary Cement Remover	<ul style="list-style-type: none"> • Efficient cleansing time (5–10) minutes • No diluting necessary—use at full strength • Color coded solution 	<ul style="list-style-type: none"> • Protects your instruments • Saves time • Economical gallon container 	<ul style="list-style-type: none"> • Removal of temporary cement including cavit and eugenol-based temporary cement • 5–10 minutes only • Place in glass beaker inside ultrasonic cleaner 	Gallon— yellow Liquid	50036212* 4000242* (Canadian customers only)
VACUKLEEN E²® EVACUATION SYSTEM CLEANER					
VacuKleen E² Evacuation System Cleaner	<ul style="list-style-type: none"> • Highly effective powdered formula • Non-corrosive and low-foaming • Biodegradable & non-toxic • Convenient, easy-mix packets • Power of protease enzymes • Fresh, clean, scent of lemon • Simple to use 	<ul style="list-style-type: none"> • Removes hard-to-clean protein-based deposits that soapy solutions can leave behind • Environmentally friendly and safe to work with • Pre-measured packets—saves time and minimizes waste • Effective and quick dissolving of blood, tissue, saliva and encrusted residue in evacuation lines • Real lemon extracts gives a clean smell that eliminates evacuation system odor and makes the operatories smell fresh • Add powder packet to warm water to make a low-foaming solution 	<ul style="list-style-type: none"> • Cleaning of your evacuation equipment • Only mix what you need so you eliminate waste and storage problems 	Box of 25 pouches Box of 60 pouches Box of 140 pouches	50036100 — Unit Dose Pouch for Quarts Intro Pack— Box of 60 Pouches (Case Lot 24) 50036101 — Unit Dose Pouch for Quarts Clinic Pack— Box of 140 Pouches (Case Lot 12) 50036102 — Unit Dose Pouch for Gallons Clinic Pack— Box of 25 Pouches (Case Lot 24)

* Each VK Solution Sold in Case Lots of Four (4) Gallons Only.